

WE RAISE THE BAR. WE SET THE STANDARD.

Project Item Number Quantity

-A-73



Ground Beef Fat Analyzer

MADE IN THE USA

Model:

> FA-73: GROUND BEEF FAT ANALYZER

Ground Beef Fat Analyzer

Features:

- Industry standard for measuring the fat content of ground beef
- Fast, simple and affordable
- Single lever manual operation
- Three minute cycle time
- Comes complete with carrying case and glass vials
- Self calibrating machine as each time you measure the fat content you are setting the machine to read zero. No external calibration is needed

Univex Corporation

Standard Accessories:

- Sample measuring device
- Carrying case
- Six glass vials

Optional Accessories:

• Carton of 48 glass vials

Warranty:

• 90 day, carry-in parts and labor




Univex Corporation www.Univexcorp.com


3 Old Rockingham Rd Salem, NH 03079 United States of America TEL: (800)258-6358 INTL TEL: (603)893-6191 FAX: (800)356-5614 INTL FAX: (603)893-1249



Ground Beef Fat Analyzer

FA-73





| Accuracy | | | | | |
|------------|-------|---------------------|--------|--|--|
| | AOAC | Modified Babcock | Univex | | |
| Regular | 23.84 | 24.7 | 24.1 | | |
| Lean | 21.58 | 21.8 | 22.1 | | |
| Extra Lean | 15.21 | 16.1 | 16.2 | | |

Accuracy: The FA-73 was tested against the AOAC method of determining fat content as well as the Modified Babcock method. The test was done with three samples of regular, lean, and extra lean meat being tested and averaged; the results are within 1%.

NEMA #: 5-15-P

| Electrical Information (60Hz) | | | | | | |
|-------------------------------|---------|-------|-------------|-----------|--|--|
| | Voltage | Phase | hp(kW) | Amps | | |
| Standard Voltage | 115 | 1 | 1.00 (0.75) | 9.0 | | |
| Optional Voltage | 220 | 1 | 1.00 (0.75) | 4.98-4.50 | | |

| Shipping Information | | | | |
|----------------------|-----------------|--|--|--|
| Net Weight | Shipping Weight | Overall Dimensions | | |
| 8 lbs. (3.6Kg) | 10 lbs. (4.5Kg) | 12" (30.5cm) H x 12" (30.5cm) W x 12" (30.5cm) D | | |



Univex Corporation www.Univexcorp.com

3 Old Rockingham Rd Salem, NH 03079 United States of America TEL: (800)258-6358 INTL TEL: (603)893-6191 FAX: (800)356-5614 INTL FAX: (603)893-1249